



Basic Tutorial

Home Module

ChromTag

Home Module

Peaks Visualization Module

Data Filtering and Grouping Module

Differential Peak Analysis Module

Differential Peak Annotation Module

Differential Peak Visualization Module

Functional Enrichment Module

Motif Enrichment Module

Help

About

Overview

Home Module

	S1	S2	S3	S4
P1	0	0	233	186
P2	0	2	166	89

Peak Count Data

Peaks Visualization

Enrichment Analysis

Differential Analysis Visualization

Motif Enrichment

Differential Peak Detection

Gene Annotation

Grouping

Volcano Plot

Heatmap

Upset Plot

Import Data

Data Source:

Example dataset

Example dataset

Upload data

(GSE145187). Each mark contains two biological replicates

START

Instruction

The **Home Module** serves as the entry point of **ChromTag**, providing an integrated interface for initializing data input and establishing the foundation for subsequent analyses. Users may choose to begin with the **system-provided example dataset** or **upload custom peak count matrices** that have been pre-merged across samples. This module also manages species selection and displays an immediate preview of the imported dataset, enabling users to verify data structure and completeness before progressing to downstream analytical modules. Together, these components ensure that the analysis pipeline begins with properly formatted and biologically consistent input data, thereby supporting accurate and reliable interpretation in later stages.

Users can choose between “**Example Dataset**” and “**Upload Data**”.

If the **Example Dataset** is selected, the application automatically loads a **pre-merged peak count matrix** containing H3K27me3 and H3K4me3 samples, each with two biological replicates.

Home Module

Import Data

Data Source:

Upload data

Select Species:

Human

Upload CSV File

Browse...

No file selected

If you choose to upload custom data, you must also select the appropriate species.

1	Chromosome	Start	End	H3K27me3_rep1	H3K27me3_rep2	H3K4me3_rep1	H3K4me3_rep2
2	chr1	27805	30656	0	0	233	186
3	chr1	135727	140083	0	2	166	89
4	chr1	198335	200968	0	0	219	169
5	chr1	391014	392146	1	3	0	1
6	chr1	392617	393317	2	1	0	0
7	chr1	393808	394451	2	1	0	0
8	chr1	440860	441530	4	2	0	1
9	chr1	442100	442830	2	0	0	0
10	chr1	492154	494761	0	1	73	36

Uploaded data should be **pre-merged across samples** and must follow a specific format. The file should be a **CSV** with a header row.

The first column must contain the chromosome name (e.g., “chr1”), the second and third columns should specify the **start** and **end** positions, and all subsequent columns should represent the **peak count values** for each sample.

Each experimental group must include **at least two biological replicates** to ensure proper statistical analysis.

Home Module

Import Data

Data Source:

Example dataset

Example dataset

Upload data

(GSE145187). Each mark contains two biological replicates , and the peaks have been pre-processed into a merged count matrix ready for analysis. You can load this dataset to explore the features and workflow of ChromTag without uploading your own files.

▶ START

After choosing the data source, users can click the **“Start”** button to load the dataset and proceed to the next module.

Data Preview

Show 10 entries

Search:

	Chromosome	Start	End	H3K27me3_rep1	H3K27me3_rep2	H3K4me3_rep1	H3K4me3_rep2
1	chr1	27805	30656	0	0	233	186
2	chr1	135727	140083	0	2	166	89
3	chr1	198335	200968	0	0	219	169
4	chr1	391014	392146	1	3	0	1
5	chr1	392617	393317	2	1	0	0
6	chr1	393808	394451	2	1	0	0
7	chr1	440860	441530	4	2	0	1
8	chr1	442100	442830	2	0	0	0
9	chr1	492154	494761	0	1	73	36
10	chr1	587572	588358	7	2	0	0

Showing 1 to 10 of 220,230 entries

Previous

1

2

3

4

5

...

22,023

Next

Once the data is successfully loaded, the system will display the **peak count matrix** in a preview table for inspection.

Peaks Visualization Module

This step is optional and can be skipped.


 Assign Sample Colors

Sample:

H3K4me3_rep2

Choose a Color:

#BE9CD6



✓ Submit


Sample	Color
H3K27me3_rep1	#6C9FB8
H3K27me3_rep2	#A587C9
H3K4me3_rep1	#7AC2B8
H3K4me3_rep2	#C5D192

Users can select a sample and assign a preferred color using the color picker. Once submitted, the chosen color will be applied consistently across the **Chromosome Coverage Plot** and the **TSS Profile Plot**, enhancing visual distinction among samples.

Peaks Visualization Module


This step is optional and can be skipped.

Users may select one or more samples and specify a chromosome to visualize peak coverage. After making the selections, click the “**Plot**” button to generate the corresponding plot.



Peaks Visualization Module

This step is optional and can be skipped.



Users can select one or more samples and define the genomic window around the transcription start site by specifying the upstream and downstream distances. After configuring the parameters, click the “**Plot**” button to generate the plot.

Data Filtering and Grouping Module

Before Filter

After Filter

Count Threshold: ?

5

✓ Submit

Show

5

entries

Search:

hromosome	Start	End	H3K27me3_rep1	H3K27me3_rep2	H3K4me3_rep1	H3K4me3_rep2
1r1	27805	30656	0	0	233	186
1r1	135727	140083	0	2	166	89
1r1	198335	200968	0	0	219	169
1r1	391014	392146	1	3	0	1
1r1	392617	393317	2	1	0	0

Showing 1 to 5 of 220,230 entries

Previous

1

2

3

4

5

...

44,046

Next

Users may specify a threshold value to filter the data.

During filtering, the system sums the read counts for each genomic region across samples and retains only those entries whose total count exceeds the defined threshold.

> 5

< 5

Data Filtering and Grouping Module

↔ Sample Grouping

Select Samples:

Select the samples that should be grouped together.

H3K27me3_rep1 × H3K27me3_rep2 ×

Group Name:

Enter a name for the current group.

H3K27me3

Add Group Clear Groups

Use the “**Add Group**” button to create the group.
Use the “**Clear Groups**” button to remove all existing groups.

Current Groups:

\$H3K27me3
[1] "H3K27me3_rep1" "H3K27me3_rep2"

All defined groups are listed under “**Current Groups.**”

✓ Submit

After all groups have been defined, click the “**Submit**” button to save the grouping information.

Differential Peak Analysis Module

↖ Differential Analysis Parameters

Direction of Change:

Both Directions ▼

Apply Independent Filtering for Low-Count Peaks

☒ Yes ☐ No

Significance Level:

0.01

P-value:

Benjamini–Hochberg (BH) ▼

▶ Run Differential Analysis

Direction of Change

This setting lets users choose whether to detect peaks that are increased or decreased between groups, allowing the analysis to focus on a specific direction of change.

Independent Filtering

If you enable independent filtering, peaks with extremely low counts will be excluded. The significance level controls the cutoff applied during this process.

Significance Level

This value specifies the statistical significance threshold used when performing independent filtering.

p-value Adjustment Method

This option allows users to select a method for multiple-testing correction.

- The Benjamini–Hochberg (BH) method controls the false discovery rate (FDR) and is well suited for analyses with many tests.
- The Holm method controls the family-wise error rate (FWER) and is more conservative, reducing false positives but potentially missing true signals.

Differential Peak Analysis Module

Comparison between two groups

Adjusted P-value Threshold:

0.05

Log2 Fold Change Threshold:

2

✓ Submit

Results can be filtered using adjusted p-values and log₂ fold-change thresholds.

Comparison between multiple groups

Differential Peak Results

1 vs 21 vs 32 vs 3

Show 10 entries

	seqnames	start	end	H3K27me3_rep1	H3K27me3_rep2
1	chr1	27805	30656	0	0
2	chr1	135727	140083	0	2

Displays multiple analysis results.

Selected Comparison:

1 vs 2

Select one comparison to use for downstream analysis.

Adjusted P-value Threshold:

0.05

Log2 Fold Change Threshold:

2

✓ Submit

Differential Peak Annotation Module

Input Annotation Parameters

Upstream (bp):

3000


Downstream (bp):

3000

▶ Run Peak Annotation

To annotate the peaks, first specify the upstream and downstream distances to define the genomic region of interest, and then click the “Run Peak Annotation” button to execute the annotation.

Differential Peak Visualization Module



Differential Peak Visualization Module

Up Genes Preview

Upregulated Genes List

DVL1

MRPL20

NOL9

DNAJC11

DHRS3

RPL11

ELOA-AS1

SRRM1

RSRP1

MACO1

Show 5 entries

Search:

	seqnames	start	end	width	strand	H3K27me3_rep1	H3K27me3_rep2	H3K4me3_
17	chr1	1346468	1350600	4133	*	0	0	
20	chr1	1405327	1408434	3108	*	0	0	
77	chr1	6552669	6555289	2621	*	0	0	
82	chr1	6699924	6702352	2429	*	0	0	
142	chr1	12614231	12620067	5837	*	0	0	

Showing 1 to 5 of 497 entries

Download

This section displays the upregulated and downregulated gene lists, as well as the peak tables generated from the filtering step.

Functional Enrichment Module

Analysis Type:

GO ▼

GO Ontology:

All ▼

p-value Cutoff:

0.05

q-value Cutoff:

0.05

Upregulated Genes List

DVL1
MRPL20
NOL9
DNAJC11
DHRS3
RPL11
ELOA-AS1
SRRM1
RSRP1
MACO1
SRM1

Downregulated Genes List


CFAP74
MEGF6
CAMTA1-AS3
CAMTA1
SLC45A1
CA6
SPSB1
LINC02606
PDPN
KAZN
KIF11

▶ Run Enrichment Analysis

Adjust the analysis type and corresponding parameters as needed.

For GO or KEGG analysis, users may enter or edit gene lists using gene symbols to define the input for enrichment testing. When GSEA is selected, the input is limited to the full set of genes annotated in the Gene Annotation module.

Functional Enrichment Module



When GSEA is selected, users must choose a specific enrichment entry before generating the corresponding GSEA plot.

Motif Enrichment Module

Upregulated Peaks List

chr3-150407290-150415877	▲
chr16-2150741-2157314	●
chr16-88972494-88978084	
chr6-119343897-119351533	
chr12-56682168-56690455	
chr8-38173859-38178882	
chr9-128688198-128692391	
chr11-65420641-65428391	
chr10-92688412-92694886	
chr6-27889944-27897392	▼
chr11-22222222-22222222	✎

Downregulated Peaks List

chr7-154746354-154754198	▲
chr11-20155235-20165969	●
chr7-27100472-27112603	
chr1-208072367-208087339	
chr1-8323134-8331215	
chr11-118150179-118161808	
chr20-47665285-47675999	
chr14-64734713-64750256	
chr1-28810214-28820062	
chr4-182814697-182830407	▼
chr12-118887818-11889817	✎

This section displays the top 200 upregulated and 200 downregulated peaks, which are selected from the differential peaks identified through the volcano plot filtering step. You can also input and modify your own lists of upregulated and downregulated peaks for further analysis.

Number of Top Transcription Factors to Plot

10

Display Motif GC Content

- ☒ Yes
- ☐ No

Enable Clustering

- ☒ Yes
- ☐ No

- Number of Top Transcription Factors to Plot:** Set how many top transcription factors you want to visualize (default is 10).
- Display Motif GC Content:** Choose whether to show the GC content of motifs.
- Enable Clustering:** Decide whether to enable clustering for grouping similar motifs in the plot.